






Zeramex XT

Zeramex Digital Solutions


















Range

XT

Zeramex Digital Solutions*	Art. No.	Name	Material
	SB15551	Customized Abutment SB (incl. Vicarbo screw)	ZrO ₂ -TZP Vicarbo
	SB15552	Customized zirconia crown for Zerabase <u>without</u> screw hole (unglazed & excl. Zerabase Abutment) SB	ZrO ₂ -TZP
	SB15553	Customized zirconia crown for Zerabase <u>with</u> screw hole (unglazed & excl. Zerabase Abutment) SB	
	SB15554	Customized zirconia Gingivaformer SB (incl. Vicarbo screw)	ZrO ₂ -TZP Vicarbo
	SB15560	Monolithic crown SB (incl. Vicarbo screw)	ZrO ₂ -TZP Vicarbo
	SB15561	2 monolithic crowns splinted SB (incl. Vicarbo screw)	
	SB15565	Monolithic bridge 3-unit SB (incl. Vicarbo screw)	
	SB15566	Monolithic bridge 3-unit SB (incl. Vicarbo screw) cutback	
	SB15590	Package customized: Zirconia Gingivaformer & Abutment SB (incl. Vicarbo screw)	
	SB15591	Package Monolithic: Zirconia Gingivaformer & crown SB (incl. Vicarbo screw)	
	RB16551	Customized Abutment RB (incl. Vicarbo screw)	










Range

	RB16552	Customized zirconia crown for Zerabase <u>without</u> screw hole (unglazed & excl. Zerabase Abutment) RB	ZrO ₂ -TZP
	RB16553	Customized zirconia crown for Zerabase <u>with</u> screw hole (unglazed & excl. Zerabase Abutment) RB	
	RB16554	Customized zirconia Gingivaformer RB (incl. Vicarbo screw)	ZrO ₂ -TZP Vicarbo
	RB16560	Monolithic crown RB (incl. Vicarbo screw)	ZrO ₂ -TZP Vicarbo
	RB16561	2 monolithic crowns splinted RB (incl. Vicarbo screw)	
	RB16565	Monolithic bridge 3-unit RB (incl. Vicarbo screw)	
	RB16566	Monolithic bridge 3-unit RB (incl. Vicarbo screw) cutback	
	RB16590	Package customized: Zirconia Gingivaformer & Abutment RB (incl. Vicarbo screw)	
	RB16591	Package Monolithic: Zirconia Gingivaformer & crown RB (incl. Vicarbo screw)	
	WB17551	Customized Abutment WB (incl. Vicarbo screw)	
	WB17552	Customized zirconia crown for Zerabase <u>without</u> screw hole (unglazed & excl. Zerabase Abutment) WB	ZrO ₂ -TZP
	WB17553	Customized zirconia crown for Zerabase <u>with</u> screw hole (unglazed & excl. Zerabase Abutment) WB	
	WB17554	Customized zirconia Gingivaformer WB (incl. Vicarbo screw)	ZrO ₂ -TZP Vicarbo
	WB17560	Monolithic crown WB (incl. Vicarbo screw)	ZrO ₂ -TZP Vicarbo
	WB17561	2 monolithic crowns splinted WB (incl. Vicarbo screw)	





Range



	WB17565	Monolithic bridge 3-unit WB (incl. Vicarbo screw)	ZrO ₂ -TZP Vicarbo
	WB17566	Monolithic bridge 3-unit WB (incl. Vicarbo screw) cutback	
	WB17590	Package customized: Zirconia Gingivaformer & Abutment WB (incl. Vicarbo screw)	
	WB17591	Package Monolithic: Zirconia Gingivaformer & crown WB (incl. Vicarbo screw)	






Zeramex Digital Solutions Multilayer*	Art. No.	Name	Material
	SB15751	Customized Abutment SB, ML (incl. Vicarbo screw)	ZrO ₂ -TZP Vicarbo
	SB15752	Customized zirconia crown for Zerabase <u>without</u> screw hole, ML (unglazed & excl. Zerabase Abutment) SB	ZrO ₂ -TZP
	SB15753	Customized zirconia crown for Zerabase <u>with</u> screw hole, ML (unglazed & excl. Zerabase Abutment) SB	
	SB15754	Customized zirconia Gingivaformer SB, ML (incl. Vicarbo screw)	ZrO ₂ -TZP Vicarbo
	SB15760	Monolithic crown SB, ML (incl. Vicarbo screw)	ZrO ₂ -TZP Vicarbo
	SB15761	2 monolithic crowns splinted SB, ML (incl. Vicarbo screw)	
	SB15765	Monolithic bridge 3-unit SB, ML (incl. Vicarbo screw)	
	SB15766	Monolithic bridge 3-unit SB, ML (incl. Vicarbo screw) cutback	
	SB15790	Package customized: Zirconia Gingivaformer & Abutment SB, ML (incl. Vicarbo screw)	



Range

	SB15791	Package Monolithic: Zirconia Gingivaformer & crown SB, ML (incl. Vicarbo screw)	ZrO ₂ -TZP Vicarbo
	RB16751	Customized Abutment RB, ML (incl. Vicarbo screw)	
	RB16752	Customized zirconia crown for Zerabase <u>without</u> screw hole, ML (unglazed & excl. Zerabase Abutment) RB	ZrO ₂ -TZP
	RB16753	Customized zirconia crown for Zerabase <u>with</u> screw hole, ML (unglazed & excl. Zerabase Abutment) RB	
	RB16754	Customized zirconia Gingivaformer RB, ML (incl. Vicarbo screw)	ZrO ₂ -TZP Vicarbo
	RB16760	Monolithic crown RB, ML (incl. Vicarbo screw)	ZrO ₂ -TZP Vicarbo
	RB16761	2 monolithic crowns splinted RB, ML (incl. Vicarbo screw)	
	RB16765	Monolithic bridge 3-unit RB, ML (incl. Vicarbo screw)	
	RB16766	Monolithic bridge 3-unit RB, ML (incl. Vicarbo screw) cutback	
	RB16790	Package customized: Zirconia Gingivaformer & Abutment RB, ML (incl. Vicarbo screw)	
	RB16791	Package Monolithic: Zirconia Gingivaformer & crown RB, ML (incl. Vicarbo screw)	
	WB17751	Customized Abutment WB, ML (incl. Vicarbo screw)	ZrO ₂ -TZP
	WB17752	Customized zirconia crown for Zerabase <u>without</u> screw hole, ML (unglazed & excl. Zerabase Abutment) WB	
	WB17753	Customized zirconia crown for Zerabase <u>with</u> screw hole, ML (unglazed & excl. Zerabase Abutment) WB	
	WB17754	Customized zirconia Gingivaformer WB, ML (incl. Vicarbo screw)	ZrO ₂ -TZP Vicarbo

Range

XT

	WB17760	Monolithic crown WB, ML (incl. Vicarbo screw)	ZrO ₂ -TZP Vicarbo
	WB17761	2 monolithic crowns splinted WB, ML (incl. Vicarbo screw)	
	WB17765	Monolithic bridge 3-unit WB, ML (incl. Vicarbo screw)	
	WB17766	Monolithic bridge 3-unit WB, ML (incl. Vicarbo screw) cutback	
	WB17790	Package customized: Zirconia Gingivaformer & Abutment WB, ML (incl. Vicarbo screw)	
	WB17791	Package Monolithic: Zirconia Gingivaformer & crown WB, ML (incl. Vicarbo screw)	

Zeramex Digital Solutions Services*	Art. No.	Name	Material
	RB16570	3D model OK/UK (without Replica) SB/RB/WB	
	RB16571	Stain & Glaze SB/RB/WB	
-	RB16580	Zeramex SB, RB & WB customised design development	-
-	RB16581	Zeramex Scan (all products)	
-	RB16582	Zeramex SB, RB & WB virtual model design based on intraoral scan data	

***Sample images/differs from original!**

Docklocs® is a registered trademark of MEDEALIS GmbH, DE.

Santoprene® is a registered trademark of Exxon Mobil Corporation, USA.

General Information

Guarantee

CeramTec Switzerland LLC offers a lifelong guarantee for implants, and a 10-year guarantee for abutments and Vicarbo screws. The Zeramex warranty conditions apply.

Certification

CE / CE₀₄₈₃ / CE₀₅₀₀ (see p. 68)

Delivery and packaging

Delivery is in accordance with the general terms and conditions (T&Cs) of CeramTec Switzerland LLC. Intact sterile packaging protects the implant from external influences and ensures sterile storage up to the printed expiration date. Zeramex XT implants and components must be stored dry in their original packaging at room temperature and protected from sunlight. Only open the packaging shortly before surgery. We recommend comprehensive clinical, radiological and statistical documentation. The inside labels (patient label) must allow traceability of the implants.

Exclusion of liability

Zeramex XT implants are part of an overall system and may be used only with the components designed

for this system. CeramTec Switzerland LLC will not be held liable for any damage arising from improper use, or from using non-original components. The general terms and conditions of CeramTec Switzerland LLC also apply.

Training

For information on courses and further education for the Zeramex XT System, please contact us at www.zeramex.com.

Material properties

All implants and abutments are made from hot-densified zirconia oxide ATZ-HIP® (HIP = Hot Isostatic Postcompaction). For reasons of quality and strength, the implants and abutments are strictly machined into their final shape from solid, hard blanks using diamond-coated tools. The workpiece does not need finishing. This allows for highly precise and reproducible production of implants and abutments with the necessary precise fit.

Zerafil implant surface

- Microstructured
- Blasted and etched
- Hydrophilic

ZrO₂ ATZ-HIP

Zirconia dioxide, ATZ (alumina-toughened zirconia) (radiopaque)

Composition:

ZrO₂ 76%, Al₂O₃ 20%, Y₂O₃ 4%

Flexural strength: >1700 MPa (on average 2000 MPa)

CTE: 9x10⁻⁶/K

ZrO₂ TZP-A

Zirconia dioxide, TZP (tetragonal zirconia polycrystal) (radiopaque)

Composition:

ZrO₂ 95%, Al₂O₃ 5%, Y₂O₃ 0.25%

Flexural strength: 1,200 MPa.

CTE: 10,5

PEEK CLASSIX

Polyether ether ketone USP Class VI (not radiopaque)

Aluminium

Aluminium (not radiopaque)

PEEK CLASSIX CW30 LSG

Short carbon fibres (CF) in a PEEK CLASSIX LSG matrix (not radiopaque)

Composition:

CF 30%, PEEK Classix LSG 70%

Flexural strength: >130 MPa

Vicarbo

Unidirectional carbon fibres (CF) in a PEEK matrix (not radiopaque)

Composition:

CF 60%, PEEK 40%

Flexural strength: >900 MPa.

General Information



	CE Mark for MDR (2017/745) Class I devices		CE Mark for legacy devices Class IIa / IIb
	CE Mark for MDR (2017/745) Class IIa / IIb devices		Caution
	Article number		Expiry date
	Batch code		Observe the instructions for use
	Do not reuse		Date of manufacture
	Non-sterile		Manufacturer
	Sterilisation by steam or dry heat		Do not use if the packaging is damaged
	Keep away from direct sunlight		
	Keep dry		Medical device
	Unique Device Identifier		EU authorised representative
	Sterile barrier system		Double packaging with single internal sterile barrier
	Do not resterilise		Gamma sterilisation
Rx only	CAUTION: US federal law restricts the sale of this product to a physician or on his or her order.		
	Implant Small ø3.5 mm SB		Implant Regular ø4.2 mm RB
	Implant Wide ø5.5 mm WB		

Contact

Europe

Order hotline for Germany, Switzerland and Austria

Phone: 00800 93 55 66 37

Fax: 00800 93 55 63 77

Email: order@zeramex.com

Zeramex Digital Solutions

Phone: 00800 04 00 13 33

Email: digitalsolutions@zeramex.com

Switzerland

CeramTec Switzerland LLC

Bodenackerstrasse 5

8957 Spreitenbach/Switzerland

Phone: 0041 44 388 36 36

Fax: 0041 44 388 36 39

Email: info@zeramex.com

Web: www.zeramex.com

Germany (EU sales)

CeramTec Dentalvertriebs GmbH

Wallbrunnstrasse 24

79539 Lörrach, Germany

Phone: +49 7621 1612 749

Fax: +49 7621 1612 780

E-mail: info@zeramex.com

Web: www.zeramex.com

Global

You can find the contact details of all sales partners on our website

www.zeramex.com/contact



Zeramex XT

For natural red-white aesthetics





