

Tightening torques ZERAMEX® screws



VICARBO® screws ZERAMEX® P6



P15450



P6 Abutment
SN Straight

Tightening torque: 25Ncm



P15451



P6 Abutment
SN Angled

Tightening torque: 25Ncm



P35550



P6 Temporary
SN

Tightening torque: 15Ncm



P16450



P6 Abutment
RN Straight

Tightening torque: 35Ncm



P16451



P6 Abutment
RN Angled

Tightening torque: 35Ncm



P36550



P6 Temporary
RN

Tightening torque: 15Ncm

Lab screws ZERAMEX® P6



P35450



P6 Abutment
SN Straight

Tightening torque: 5Ncm



P35451



P6 Abutment
SN Angled

Tightening torque: 5Ncm



P36450



P6 Abutment
RN Straight

Tightening torque: 5Ncm



P36451



P6 Abutment
RN Angled

Tightening torque: 5Ncm